

Whole Genome Bisulfite Sequencing Protocol

This protocol combines NEBNext oligos, the KAPA Hyper Prep kit, and a bisulfite conversion kit to generate WGBS libraries. Please look at the manufacturer's instructions, product information, and referenced papers for more details.

The following kits and reagents are required:

- i) QIAquick PCR Purification Kit (QIAGEN, cat# 28104)
- ii) NEBNext® Multiplex Oligos for Illumina® (Methylated Adaptor, Index Primers Set 1) (NEB cat# E7535S)
- iii) KAPA Hyper Prep Kit (ROCHE, cat# KK8500, KK8502, KK8504)
- iv) Bisulfite conversion kit (ex. Zymo Gold)
 - ❖ There is no specific conversion kit recommended due to differing conversion rates seen between experiments. Optimization is recommended.
- v) PfuTurbo Cx Hotstart DNA Polymerase (Aligent, cat# 600410)
- vi) KAPA Library Quantification Kit (ROCHE, cat# KK4854)
- vii) AMPure XP beads (Beckman, cat# A63881)
- viii) Fragment Analyzer kit DNF-474-0500
- ix) Qubit™ dsDNA HS (ThermoFisher, cat# Q32854)

I. Preparation genomic DNA

Use a gDNA extraction protocol comparable for bisulfite conversion like Bradley's or a commercial kit like QIAGEN.

II. Fragmentation

Optimum size of the fragmented gDNA is around 250bp (peak size). Use a Covaris S220/E220 sonicator to shear DNA into small fragments.

- 1) Prepare 5 µg of gDNA.
- 2) Mix the gDNA samples with 0.5% W/W (25 ng) lambda DNA and top up to a final volume of 130 µl in TE buffer or 10mM Tris-Cl (pH8.0).
- 3) Fragment the DNA using the following settings:

Peak Incident Power	175W
Duty Factor	10%
Cycle per burst	200
Time	3mins, 30 seconds on 30 seconds off
Water level	12.5
Temperature	4°C – 6°C
Volume	130 µl

- 4) Verify the size of the fragmented DNA on a 1.2% gel or on a Fragment Analyzer using the DNF-474 kit.
- 5) Purify the DNA by using QIAquick PCR purification Kit and elute in 52µl of EB or 10mM Tris-Cl (pH8.0).

III. End Repair and A-Tailing

The following uses the KAPA Hyperprep kit (Cat# KK8502).

- 6) Prepare the following reaction mix on ice:

End Repair & A-Tailing Mix	
Reagent	Volume (µL)
Fragmented DNA Sample (1ug)	50
End Repair & A-Tailing Buffer	7
End Repair & A-Tailing Enzyme Mix	3
Total Volume	60

- 7) Vortex gently and spin down briefly.
- 8) Incubate the reaction in a thermocycler for 30 minutes at 60 °C, lid at 85 °C, store at 4 °C.
- 9) Proceed to Adapter ligation.

IV. Methylated Adapter Ligation

The following reaction volume works for 100 ng – 1 µg of DNA.

Methylated adapter MUST be used to avoid conversion of “C” in the adapter to “T” during bisulfite conversion in the following step. CT conversions occurring within the adapter will lead to a change in the adapter sequence resulting in a failure of the final library sequencing.

- ❖ Bring AMPure XP beads to room temperature before using.

10) Prepare the following reaction mix on ice:

Ligation Mix	
Reagent	Volume (µL)
DNA Sample	60
Ligation buffer	30
Methylated Adapter from NEB	2.5
DNA Ligase	10
PCR-grade Water	7.5
Total Volume	110

- 11) Incubate the reaction for 15 minutes at 20°C.
- 12) Add 3 µl of USER enzyme to the ligation mixture and mix well.
- 13) Incubate the reaction for 15 minutes at 37 °C.
- 14) Purify the DNA with 1x AMPure XP beads (113 µl).

V. Bisulfite Conversion of DNA

- 15) Follow the procedure of corresponding Bisulfite conversion kit (such as Zymo Gold).
- 16) Elute the converted DNA in 10 µl of elution buffer.
 - ❖ SAFE STOPPING POINT – Store DNA in -80°C.

VI. Amplification of BS converted library

17) Prepare the following reaction:

PCR Mix	
Reagent	Volume (µL)
BS converted DNA	10µl
10x Pfu CX Buffer	5µl
10uM NEBNext Universal PCR Primer	2.5µl
10uM NEBNext Index PCR Primer	2.5µl
10mM dNTP	1.5µl
Milli-Q / ddH2O	27.5µl
Pfu Turbo CX Polymerase	1µl
Total Volume	50µl

18) Set up the following cycling condition on the thermal cycler:

Thermocycler Conditions			
Stage 1	95°C	2 minutes	
	98°C	30 seconds	
Stage 2	98°C	15 seconds	Repeat 4-8 cycles
	60°C	30 seconds	
	72°C	4 minutes	
Stage 3	72°C	10 minutes	
	4°C	Hold	

- ❖ 4-8 cycles should use for the amplification. If the DNA concentration use for bisulfite conversion is low, more cycles should be use.

19) Double size select DNA with 0.7x AMPure XP beads (keeping the eluent), followed by 0.2x purification.

20) Elute in 20 µL EB

21) Quantify the DNA using both Qubit and KAPA Library Quantification Kit.

22) Run 5ng of DNA on the Fragment Analyzer using DNF-474 analysis to check the average size of library prepared.

23) If necessary, pool libraries together for sequencing.

- ❖ A considerable amount of ssDNA might be present in the final library due to the DNA damage during bisulfite conversion. Quantification of the library should be performed cautiously, taking into account the existence of ssDNA.