

μChIP

This protocol is modified from John Arne Dahl and Inkyung Jung, *et al.*, 2016 to prepare ChIP samples from low cell numbers. Certain parts of the protocol need to be optimized depending on the antibody and cell types used. Please look at the original paper for further details.

Please look at the separate protocol for library preparation if continuing to Illumina sequencing

Crosslinking buffer				
Reagents	Stock concentration	Final concentration	5 ml	10 ml
Formaldehyde	16%	1.1%	343.75 μl	687.5 μl
PIC (protease inhibitor)	25X	1X	200 μl	400 μl
Sodium butyrate	0.5M	20 mM	200 μl	400 μl
PMSF	0.2M	1 mM	25 μl	50 μl
EDTA	0.5M	1 mM	10 μl	20 μl
PBS	---	---	4221.25 μl	8442.5 μl

NOTES:

- ❖ freshly prepared, store at RT
- ❖ formaldehyde waste should be disposed in separate container
- ❖ PIC, Sodium butyrate and PMSF should be added last
- ❖ PIC, Sodium butyrate and PMSF has 1 month shelf life
- ❖ Formaldehyde should be methanol free
- ❖ Dispose sodium butyrate once opened
- ❖ PIC pellet in 2mL is 25X

Crosslinking wash buffer				
Reagents	Stock concentration	Final concentration	5 ml	10 ml
Sodium butyrate	0.5M	20 mM	200 μl	400 μl
PIC	25X	1X	200 μl	400 μl
PMSF	0.2M	1 mM	25 μl	50 μl
EDTA	0.5M	1 mM	10 μl	20 μl
PBS	---	1X	4565 μl	9130 μl

- ❖ Prepare buffer right before use. Store at 4C

RIPA buffer				
Reagents	Stock concentration	Final concentration	5 ml	10 ml
Triton X-100	10%	1%	500 μl	1000 μl
NaCl	5M	140 mM	140 μl	280 μl
Tris-HCl (pH 8.0)	1M	10 mM	50 μl	100 μl
SDS	10%	0.1%	50 μl	100 μl
Na-deoxycholate (DOC)	10%	0.1%	50 μl	100 μl
EDTA	0.5M	1 mM	10 μl	20 μl
EGTA	0.5M	0.5 mM	5 μl	10 μl
ddH ₂ O (Milli-Q)	---	---	4195 μl	8390 μl

Lysis buffer				
Reagents	Stock concentration	Final concentration	1ml	5ml
SDS	10%	0.8%	80 μl	400 μl
Na butyrate	0.5M	20mM	40 μl	200 μl
PMSF	0.2M	1mM	5 μl	25 μl
PIC	25X	1X	40 μl	200 μl
PBS	---	---	835 μl	4175 μl

RIPA dilution buffer				
Reagents	Stock concentration	Final concentration	5 ml	10 ml
Triton X-100	10%	1.25%	625 μl	1250 μl
Sodium butyrate	0.5M	20 mM	200 μl	400 μl
PIC	25X	1 X	200 μl	400 μl
NaCl	5M	175 mM	175 μl	350 μl
Na-deoxycholate (DOC)	10%	0.125%	62.5 μl	125 μl
Tris-HCl (pH 8.0)	1M	10 mM	50 μl	100 μl
PMSF	0.2M	1 mM	25 μl	50 μl
EDTA	0.5M	1 mM	10 μl	20 μl
EGTA	0.5M	0.625 mM	6.25 μl	12.5 μl
ddH ₂ O (Milli-Q)	---	---	3646.25 μl	7292.5 μl

ChIP elution buffer				
Reagents	Stock concentration	Final concentration	5 ml	10 ml
SDS	10%	1%	500 μl	1000 μl
Tris-HCl (pH 7.5)	1M	20 mM	100 μl	200 μl
RNaseA	10 mg/mL	30 μg (per 150 μl)	100 μl	200 μl
NaCl	5M	50 mM	50 μl	100 μl
EDTA	0.5M	5 mM	50 μl	100 μl
ddH ₂ O (Milli-Q)	---	---	4200 μl	8400 μl

Note: Add RNaseA separately to individual tube and not to master mix.

I. Crosslinking of cells

- 1) Add 50 μl of PBS with 2% formaldehyde to 50 μl of cells to get a 1% final formaldehyde concentration and vortex (gently!), incubate at RT for 8 min, and vortex once more.
 - ❖ mESCs were prepared at 1 million cells per 50 μl
- 2) Add 5.3 μl of 2.5 M glycine stock (final conc. 125 mM), mix by gently vortexing, incubate at RT for 5 mins, and vortex once during the incubation.
- 3) Centrifuge the crosslinked cells at (centrifuge speed is dependent on cell type) at 4°C in a swing-bucket rotor with soft deceleration.
- 4) Wash the pellet twice with 400 μl ice cold PBS. Snap freeze after the last wash, store at -80°C.

II. Bead preparation

- 5) Transfer 2 μl of Dynabeads Protein A (per IP) to a 1.5 ml tube.
- 6) Place on a magnetic rack.
- 7) Remove excess buffer and wash with 500 μl RIPA buffer.

- 8) Resuspend beads in RIPA buffer to 100 μl per IP, move to new 200 μl tube.
- 9) Add 2 μg antibody to the tube.
- 10) Incubate at 40 rpm on rotator for at least 4 hrs at 4°C.

III. Chromatin preparation

- 11) Slowly thaw a tube(s) of frozen cell pellet to ice.
- 12) Add 100 μl of crosslinking buffer to the pellet.
- 13) Incubate at room temperature for 10 mins and vortex gently twice.
- 14) Add 3.5 μl of 2.5 M Glycine, vortex gently and incubate for 5mins.
- 15) Centrifuge for (centrifuge speed differs for cell type) at 4°C.
- 16) Wash the pellet twice with 400 μl of crosslinking wash buffer.
- 17) Discard the supernatant and resuspend in 120 μl lysis buffer.
- 18) Transfer to sonication tube.
- 19) Sonicate at optimized conditions.

❖ For example, Covaris S220/E220 for tube of 1 million mESC

Peak Incident Power	175.0
Duty Factor	10.0
Cycle per burst	200
Time (seconds)	400
Water level	12.5
Temperature	2°C – 7°C
Volume	130 μl

- 20) Transfer lysate to new 1.5 ml tube, rinse sonication tube with 170 μl RIPA dilution buffer.
- 21) Centrifuge for 12,000 x g for 10 mins at 4°C.
- 22) Transfer supernatant to a new 1.5 ml tube.

IV. IP

- 23) Quick spin the pre-incubated beads and antibody complex to collect liquid at the cap.
- 24) Pellet on cold magnetic rack (keep on ice).
- 25) Remove the supernatant.
- 26) Wash in 130 μl RIPA buffer, vortex roughly.
- 27) Repeat the washing for a total of two washes.
- 28) Transfer appropriate volume of chromatin to each IP tube, keep the same volume as input. Also keep a small amount for size check.
 - ❖ Volume per of chromatin depends on how many marks you are doing and how you are splitting a tube of cells (ie. 1 mark and 1 input, roughly 65 μl)
- 29) Add 4 μg crosslinked recombinant histone octamers and 1.25 μg of non-immunized rabbit IgG to each ChIP reaction.
- 30) Top up IP volume to ~180 μl with RIPA dilution buffer.
- 31) Rotate the tubes at 4°C for at least 30hrs.

V. Washing

- 32) Wash the chromatin-antibody-bead complex 4 times with 100 μl cold RIPA buffer.
 - ❖ Concentrations of SDS and NaCl can be optimized for the antibody.
- 33) Wash the beads in 100 μl TE buffer and transfer to a new 1.5 ml tube.
 - ❖ Chill the TE before use

VI. Elution

- 34) Pellet on cold magnetic rack, remove the TE buffer and add 150 μl CHIP elution buffer.
- 35) Incubate at 37°C for 1hr, 1000 rpm on thermomixer.
- 36) Add 15 μl of Proteinase K (20 mg/ml) and incubate at 60°C for 1 hrs, 1000rpm, followed by 68 °C for 3hrs
- 37) Transfer eluate to a new 1.5 ml tube.
- 38) Add 150 μl fresh CHIP elution buffer to the beads and further elute for 5mins.
- 39) Pool the eluate from both elution.
- 40) Purify the DNA by phenol-chloroform isoamylalcohol extraction (ethanol precipitated with 10 μl acrylamide carrier) or use a kit such as QIAGEN MinElute PCR purification kit.
- 41) Dissolve or elute in 20 μl of EB (or 10 mM Tris-HCl).

References:

Dahl, J., Jung, I., Aanes, H. et al. Broad histone H3K4me3 domains in mouse oocytes modulate maternal-to-zygotic transition. *Nature* **537**, 548–552 (2016). <https://doi.org/10.1038/nature19360>